

EDGEWOOD CHEMICAL BIOLOGICAL CENTER

U.S. ARMY RESEARCH, DEVELOPMENT AND ENGINEERING COMMAND
Aberdeen Proving Ground, MD 21010-5424

ECBC-TR-1465

GC-MS/MS ANALYSES OF BIOLOGICAL SAMPLES IN SUPPORT OF EVALUATION OF TOXICITY ASSOCIATED WITH INTRAVENOUS EXPOSURE TO VX STEREOISOMERS IN GUINEA PIGS

Jeffrey M. McGuire Linnzi K. Wright Robert L. Kristovich

RESEARCH AND TECHNOLOGY DIRECTORATE

Michael W. Busch

EXCET, INC. Springfield, VA 22151-2110

July 2017

Approved for public release: distribution unlimited.



Disclaimer	
The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorizing document	ts.

	REPORT [OCUMEN.	TATION PAGE	<u> </u>	Form Approved OMB No. 0704-0188
Public reporting burden of needed, and completing burden to Department of Respondents should be	for this collection of inform and reviewing this collect Defense, Washington He aware that notwithstandin	nation is estimated to avera ion of information. Send of eadquarters Services, Dire	age 1 h per response, including the comments regarding this burden es ctorate for Information Operations w, no person shall be subject to ar	e time for reviewing instruction stimate or any other aspect of and Reports (0704-0188), 1	ons, searching existing data sources, gathering and maintaining the data fithis collection of information, including suggestions for reducing this 215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302. It with a collection of information if it does not display a currently valid
1. REPORT DATE		2. REPORT T			3. DATES COVERED (From - To)
XX-07-2017		Final			Jan 2016 – Aug 2016
4. TITLE AND SU		ogical Samples	in Support of Evalua	ntion	5a. CONTRACT NUMBER
of Toxicity As	sociated with Ir	•	osure to VX Stereoise		5b. GRANT NUMBER
in Guinea Pigs					5c. PROGRAM ELEMENT NUMBER
6. AUTHOR(S)					5d. PROJECT NUMBER
• •	ev M · Wright	Linnzi K · Kris	tovich, Robert L. (EC	CBC): and	CB10168
Busch, Michae	•	Limizi IX., IXIIS	tovicii, Robert E. (Ex	CDC), and	5e. TASK NUMBER
,	, ,				EC WORK UNIT NUMBER
					5f. WORK UNIT NUMBER
		NAME(S) AND ADD			8. PERFORMING ORGANIZATION REPORT NUMBER
			G, MD 21010-5424 Springfield, VA 2215	51-2110	ECBC-TR-1465
			ND ADDRESS(ES)	2110	10. SPONSOR/MONITOR'S ACRONYM(S)
Defense Threa	t Reduction Ag		J. Kingman Road, N	MSC 6201,	DTRA
Fort Belvoir, V	/A 22060-6201				11. SPONSOR/MONITOR'S REPORT NUMBER(S)
12. DISTRIBUTIO	N / AVAILABILITY	STATEMENT			
	•	listribution unli	mited.		
13. SUPPLEMEN	TARY NOTES				
14. ABSTRACT:					
					metry analyses of blood, tissues,
					gs intravenously exposed to various
	•	1 1 2			X) stereoisomers. The amount of
		· / I	present in these biolo	gical samples wa	s quantified using chemical ionization
mass spectrom	etry with isotop	e dilution.			
15. SUBJECT TE	-				
			etry (GC–MS/MS)		oisomers
		oridate (VX-G)			pe dilution
<i>O</i> -ethyl <i>S</i> -(2-	diisopropylamii	noethyl) methyl	phosphonothioate (V	(X) Intra	venous exposure
16. SECURITY CL	ASSIFICATION O	F:	17. LIMITATION OF	18. NUMBER OF	19a. NAME OF RESPONSIBLE PERSON
			ABSTRACT	PAGES	Renu B. Rastogi
a. REPORT	b. ABSTRACT	c. THIS PAGE	7 77 7	60	19b. TELEPHONE NUMBER (include area code)
U	U	U	UU	60	(410) 436-7545

(410) 436-7545 Standard Form 298 (Rev. 8-98) Prescribed by ANSI Std. Z39.18 Blank

PREFACE

The work described in this report was authorized under project no. CB10168 and Institutional Animal Care and Use Committee protocol 16-472. The work was started in January 2016 and completed in August 2016, as recorded in U.S. Army Edgewood Chemical Biological Center (ECBC; Aberdeen Proving Ground, MD) Notebook 14-0084.

The use of either trade or manufacturers' names in this report does not constitute an official endorsement of any commercial products. This report may not be cited for purposes of advertisement.

This report has been approved for public release.

Acknowledgments

The authors acknowledge Ruth W. Moretz, Jeffry S. Forster, Julie A. Renner, and Bernardita I. Gaviola (ECBC) for their efforts in working with the exposures and sample collections.

Blank

CONTENTS

1.	INTRODUCTION	
2		
2.	METHODS	1
2.1	Animal Exposures	1
2.1.1	Animals	1
2.1.2	Range-Finding and LD ₅₀ Studies	2
2.1.3	Pharmacokinetic Studies	
2.2	Sample Preparation and Analysis	
2.2.1	Chemical Materials	
2.2.2	Stock Solutions and Calibration Standards	
2.2.3	Analytical Method	3
2.2.4	Sample Preparation	3
3.	RESULTS AND DISCUSSION	4
3.1	Range-Finding Studies	4
3.2	LD ₅₀ Studies	
3.3	Pharmacokinetic Studies	
	LITERATURE CITED	47
	ACRONYMS AND ABBREVIATIONS	49

FIGURES

1.	The S _N 1 reaction: nucleophilic reagent attacks both back and	
	front sides, resulting in racemization	6
2.	VX-G concentration in blood as a function of the dose	
	of individual VX stereoisomers or the racemic mixture	16
3.	VX-G concentration in blood as a function of time	
	for three different doses of P(+)-VX	45
	TABLES	
1.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(+)-VX for Range-Finding Studies	5
2.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(–)-VX for Range-Finding Studies	5
3.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(+)-VX for LD ₅₀ Studies	7
4.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(–)-VX for LD ₅₀ Studies	11
5.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to Racemic VX for LD ₅₀ Studies	14
6.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies	18
7.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies	26
8.	Results from VX-G Assays of Blood, Tissues, and Organs Following	
	Intravenous Exposure to Racemic VX for Pharmacokinetic Studies	34
9.	Pharmacokinetic Parameters from Serial Samplings Following	
	Intravenous Exposure to Various Doses of P(+)-VX, P(-)-VX, and a	
	Racemic Mixture of VX	45

GC-MS/MS ANALYSES OF BIOLOGICAL SAMPLES IN SUPPORT OF EVALUATION OF TOXICITY ASSOCIATED WITH INTRAVENOUS EXPOSURE TO VX STEREOISOMERS IN GUINEA PIGS

1. INTRODUCTION

Many organophosphorus chemical warfare nerve agents have an asymmetric phosphorous atom, such that the synthesis of these agents results in at least two stereoisomers, designated P(+) and P(-). Although in vitro studies have demonstrated that the P(-)-isomers of these agents are several orders of magnitude more potent inhibitors of acetylcholinesterase than are the P(+)-isomers (Nordgren et al., 1984; Benschop and de Jong, 1988; Ordentlich et al., 2004), only a few in vivo studies have been conducted with optically pure stereoisomers (reviewed in Benschop and de Jong, 1988, 2001). By subcutaneously exposing mice to each of the stereoisomers of sarin, soman, tabun, and O-ethyl S-(2-diisopropylaminoethyl) methylphosphonothioate (VX), it was determined that the median lethal dose (LD50) values of the P(-)-isomers of these agents were approximately half that of the racemic mixture. However, mice are not the ideal animal model for studying the toxicity associated with nerve agent exposure because they have relatively high levels of carboxylesterase activity (Maxwell et al., 1987).

The Evaluation of the Toxicity Associated with Intravenous Exposure to the Stereoisomers of Chemical Warfare Nerve Agents in Guinea Pigs, which is Institutional Animal Care and Use Committee (IACUC) protocol number 16-472 (IACUC, 2015), provides data to determine whether the individual stereoisomers of chemical warfare agents have different toxic responses in guinea pigs and to characterize the pharmacokinetics of the individual stereoisomers and their racemic mixtures. This report details the results of gas chromatography—tandem mass spectrometry (GC–MS/MS) analyses of blood, tissues, and organs that were performed to quantify the amounts of *O*-ethyl methylphosphonofluoridate (VX-G) present in guinea pigs after intravenous exposure to P(+)-VX, P(-)-VX, and a racemic mixture of VX.

2. METHODS

2.1 Animal Exposures

2.1.1 Animals

Adult, male guinea pigs weighing 350–400 g that had been surgically implanted with jugular vein catheters were purchased from Charles River Laboratories International, Inc. (Kingston, NY). Guinea pigs were single-housed in temperature- and humidity-controlled rooms (21 \pm 1 °C and 30–70%, respectively). Lights were turned on at 0600 and off at 1800. Food and water were provided ad libitum, and guinea pigs had access to enrichment items such as huts and chew toys.

2.1.2 Range-Finding and LD₅₀ Studies

Optically pure stereoisomers were separated from the racemic mixture of VX by workers at the Agent Chemistry Branch (U.S. Army Edgewood Chemical Biological Center [ECBC]; Aberdeen Proving Ground, MD) using methods described by Bae and Winemiller (2016). Guinea pigs were intravenously exposed (through their catheters) to one of the stereoisomers or the racemic mixture. A range-finding study (n = 4 or 5 animals per agent) was conducted with P(+)- and P(-)-VX to generate starting doses for the LD₅₀ studies, whereas the LD₅₀ value reported by Shih and McDonough (2000) was used as the starting dose for racemic VX. Initially, 2-propanol (IPA) was used as the solvent; however, this was changed to saline midway through the LD₅₀ studies (exposure dates on and after 10 March 2016). In the LD₅₀ studies for which IPA was the solvent (n = 9 or 10 animals per dose; injection volume of 0.5 mL/kg), the doses ranged from 100 to 1000 μ g/kg for P(+)-VX, 1.8 to 3.5 μ g/kg for P(-)-VX, and 7.0 μg/kg for racemic VX. In the LD₅₀ studies for which saline was the solvent (n = 4-10 animals per dose; injection volume was 0.5 mL/kg), the doses ranged from 175 to280 μ g/kg for P(+)-VX, 3.0 to 5.0 μ g/kg for P(-)-VX, and 3.7 to 7.0 μ g/kg for racemic VX. Toxic signs were recorded continuously for 2 h post-exposure and then intermittently until the close of business. At 24 h post-exposure, guinea pigs were euthanized with a barbiturate solution containing 100 mg/kg of sodium pentobarbital. Biosamples (blood, brain, heart, liver, lung, kidneys, and urine) were collected at time of death or euthanasia.

2.1.3 Pharmacokinetic Studies

A subset of guinea pigs (n = 4 animals per dose) that were surgically implanted with double jugular vein catheters were intravenously exposed via their left catheters to the same doses of VX that were used in the LD₅₀ studies with the saline solvent. Toxic signs were recorded as for range-finding studies, and guinea pigs were euthanized at 24 h post-exposure. Blood samples were collected via the right catheters at 0, 1, 10, 20, 30, 40, 50, and 60 min and 3, 6, and 24 h, and biosamples were collected at time of death or euthanasia.

2.2 Sample Preparation and Analysis

2.2.1 Chemical Materials

EA 1207 (VX-G) and deuterated (2H_5) EA 1207 were obtained from ECBC stock. Before use, the EA 1207 was verified by quantitative ^{31}P NMR spectrometry to be 78.43 ± 0.56 wt % (ECBC laboratory notebook [NB] 11-0003-114), and the 2H_5 EA 1207 was verified to be 68.81 ± 0.9 wt % (NB 11-0003-115). Potassium fluoride (KF), IPA, ethyl acetate, glacial acetic acid, and anhydrous sodium sulfate were obtained at $\geq 99\%$ purity from Sigma-Aldrich (St. Louis, MO). Sodium acetate at >99% purity was purchased from Fischer Chemicals (Fair Lawn, NJ). Ammonia and methane were obtained from Sigma-Aldrich, and helium was obtained from Messer (Malvern, PA); all were at >99.9% purity.

2.2.2 Stock Solutions and Calibration Standards

Stock solutions of VX-G and 2H_5 VX-G (internal standard [IS]) were prepared in IPA at concentrations of 1.524 mg/mL (NB 11-0003-110-01) and 1.473 mg/mL (NB 11-0003-111-01), respectively, and stored at -20 °C until use. Working solutions (5–10 μ g/mL) were prepared by diluting the stock solutions in ethyl acetate. Calibration standards of VX-G were prepared by diluting the working solutions to obtain the following 12 concentration points: 0.5, 1, 5, 10, 25, 50, 100, 200, 400, 600, 800, and 1000 ng/mL (NBs 11-0003-117-04 through 11-0003-117-15). Each calibration standard also contained 200 ng/mL of 2H_5 VX-G diluted from the working solution. All calibration standards were stored at -20 °C until analysis.

For the VX-G assays, calibration curves were constructed using the 12 calibration standards from the respective analytes, where relative response (defined as $area_{analyte}/area_{IS}$) was plotted against relative concentration (defined in nanograms per milliliter as concentration_{analyte}/concentration_{IS}). For the VX-G calibration curve, a quadratic curve fit was used with a 1/x weighting factor. Typically, these calibration curves yield correlations of $R^2 = 0.999$ over 3 orders of magnitude, where R^2 is defined as the coefficient of determination.

2.2.3 Analytical Method

VX-G sample assays were performed using an Agilent 7000A Triple Quad GC/MS instrument (Agilent Technologies; Santa Clara, CA). Gas chromatographic separations were achieved using an RTx-1701 column (30 m \times 0.25 mm i.d., 0–25 µm film thickness; Restek Corporation; Bellefonte, PA). The carrier gas was helium, and the flow rate was 1 mL/min. Injections of 2.0 or 3.0 µL were made using an Agilent 7693 ALS autoinjector into a splitless injector port at a temperature of 225 °C. The initial oven temperature of 35 °C was held for 6 s, then ramped to 100 °C at 15 °C/min, and ramped again at 35 °C/min to 175 °C. After each analysis was complete, the column was back-flushed at 280 °C for 4 min at reduced inlet pressure (–6.3 mL/min). The typical retention time for VX-G and its deuterated standard is 5.5 min.

Samples were ionized by positive-ion chemical ionization (CI) with ammonia reagent gas. CI source conditions were optimized using Fluoroether E3 tuning compound (Chemical Abstracts Service [CAS] registry number 3330-16-3; Agilent Technologies) with methane reagent gas. Mass spectra were obtained at a dwell time of 0.2 s for each transition in the multiple reaction monitoring (MRM) mode. Helium was used as the collision gas, and the collision energy (CE) was 12 V. The CE was optimized for the mass-to-charge ratio (m/z) 144 > 99 transition for VX-G and the m/z 149 > 100 transition for 2H_5 VX-G. The MassHunter software provided with the Agilent 7000A system was used to process and analyze the data. The software provides automated peak detection, calibration, and quantitation.

2.2.4 Sample Preparation

Sample preparations for this study were similar to those previously described by McGuire et al. (2015). Upon arrival at ECBC, all biological samples were stored at -80 °C until analysis. Whole blood samples were extracted for VX-G using Oasis HLB 30 μ m solid-phase

extraction (SPE) cartridges (Waters Corporation; Milford, MA), which were first conditioned with 1 mL each of ethyl acetate, IPA, and pH 4.0 acetate buffer (0.01 M sodium acetate and 0.2 M glacial acetic acid). A sample of blood in a 2.0 mL microcentrifuge tube (Sigma-Aldrich) was weighed, and then 1 mL of acetate buffer, 200 μ L of KF solution (6 M), and 1 μ L of IS, 2 H $_5$ VX-G were added. The mixture was vortexed for 10–20 s and centrifuged at 15,000 rpm for 5 min using a Micromax microcentrifuge (Thermo IEC; Needham Heights, MA). The supernatant liquid was transferred to the SPE cartridge, and the sediment at the bottom of the microcentrifuge tube was resuspended with 750 μ L of acetate buffer and 200 μ L of KF solution. This mixture was also vortex-mixed and centrifuged, and the resulting liquid was added to the original solution. After the mixture was added to the SPE cartridge, it was allowed to drain under a gentle vacuum. The analytes were eluted with 1 mL of ethyl acetate, which was collected and dried over anhydrous sodium sulfate. The ethyl acetate was withdrawn from the collection tube, filtered through a 0.2 μ m nylon Acrodisc syringe filter (Pall Gelman Laboratory; Ann Arbor, MI) into a GC autosampler vial (Agilent Technologies), and then concentrated to 50 μ L for analysis.

Tissue and organ sample extracts were prepared in a similar manner; freeze-fracture pulverization under cryogenic temperatures was performed before SPE extraction. A CryoPrep system (Covaris; Woburn, MA) was used to pulverize 0.5–1 g of tissue. The pulverized sample was mixed with 1 mL of acetate buffer, 200 μL of KF solution, and 1 μL of IS. This sample was then subjected to focused acoustics using an S-series focused ultrasonicator. This system directs precisely controlled cavitation and acoustic streaming to a focal point within a sample-treatment vessel in a noncontact, isothermal process. After centrifugation at 4500 rpm for 15 min using a Sorvall Legend X1R centrifuge (Thermo Fisher Scientific; Waltham, MA), the supernatant liquid was transferred to the SPE cartridge, and the sediment at the bottom of the sample tube was resuspended with 750 μL of acetate buffer and 200 μL of KF solution. This mixture was vortex-mixed and centrifuged, and the resulting liquid was added to the original solution. Additional sample processing was performed in a manner identical to that used for the blood samples.

3. RESULTS AND DISCUSSION

3.1 Range-Finding Studies

The following results were recorded in ECBC NB 14-0084. Table 1 summarizes the data from VX-G assays of whole blood, tissues, and organs that were obtained from guinea pigs after intravenous exposure to various doses of P(+)-VX. Table 2 shows similar results from intravenous exposures to various doses of P(-)-VX. The amounts of VX-G reported are representative of both free VX, which may be present, and VX bound to proteins and subsequently released as VX-G.

Although the guinea pigs were exposed to optically pure stereoisomers of VX, a racemic mixture of VX-G was obtained from all of the samples analyzed. Upon addition of KF to the acetate buffer solution, the nucleophilic substitution of F^- proceeds via an S_N1 reaction in which all stereochemistry is lost, and racemization occurs as shown in Figure 1.

5

Table 1. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Range-Finding Studies

Guinea	Exposure Date	Dose (mg/kg)	Sample Time	VX-G (ng/g)						
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
2		10.0*	5	5,847.3965	13,321.9823	21,393.6667	16,033.5775	22,312.9662	40,605.6615	
8		1.0*	1,440	3.4910	3.2014	11.1297	106.2745	188.5999	4.4013	
9	27 Jan 2016	5.0*	6	3,105.2312	4,543.2631	17,736.0211	4,986.2676	11,475.2847	8,114.1826	
10		2.5*	1	405.8055	1,486.9555	3,778.3479	384.9394	67.7238	445.9136	
32		0.5*	2	369.4611	4,287.1277	9,959.7540	657.6626	41.9359	441.4865	
215		0.28	1,440	4.8094	16.8223	57.0460	280.9087	933.2821	13.3417	
222	12 Apr 2016	0.56	26	148.9941	166.1160	245.4801	1,503.0742	2,960.5377	141.3365	
229	12 Apr 2016	0.23	1,440	6.4527	10.8781	67.2446	92.3158	318.3802	7.6369	
232		0.14	1,440	2.7106	1.4469	6.6270	32.1262	93.3448	5.0333	

^{*2-}Propanol was used as the solvent for these dosings.

Table 2. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Range-Finding Studies

Guinea	Exposure	Dose*	Sample Time		VX-G (ng/g)						
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
80		10.0	18	10.9380	8.7204	18.5904	4.2004	18.3587	13.8842		
110	1 Mar 2016	1.0	1,440	0.6923	0.3571	0.5751	0.4312	BDL	1.6282		
118	1 Mar 2016	3.5	93	3.3024	1.8911	3.9071	10.2575	12.3515	8.7529		
126		2.2	1,440	1.0952	0.4913	0.8004	0.6840	BDL	3.6986		

^{*2-}Propanol was used as the solvent for all dosings.

BDL, below detection limit (<0.05 ng/g).




Figure 1. The S_N1 reaction: nucleophilic reagent attacks both back and front sides, resulting in racemization.

3.2 LD₅₀ Studies

Table 3 summarizes the data from the VX-G assays of whole blood, tissues, and organs that were obtained from guinea pigs after intravenous exposure to various doses of P(+)-VX. Tables 4 and 5 show similar results from intravenous exposures to various doses of P(-)-VX and a racemic mixture of VX, respectively.

Figure 2 is a plot of the VX-G concentration in the blood as a function of the dose of the individual VX stereoisomers or the racemic mixture. Each point represents an average (n=4-10 animals) of the individual values determined at the time of death or at 24 h for each dose prepared in saline. Error bars represent the standard deviation of each average determination. Curves were fitted using a Hill-type, three-parameter pharmacodynamic model of the form $y = a \times x^b/(c^b + x^b)$, where x and y represent the corresponding axis values on Figure 2, a is the maximum VX-G concentration observed, b is a sigmoidicity factor representative of the curve shape, and c is the concentration or dose that yields a 50% effect. These dose-response curves clearly show the increased toxicity of the P(-)-VX stereoisomer and the racemic mixture versus the P(+)-VX stereoisomer. From visual inspection of the curves, the dose at which a 50% response (VX-G concentration) occurred is estimated as 4–5 μ g/kg for the P(-)-VX stereoisomer, 5–6 μ g/kg for the racemic mixture, and 260–280 μ g/kg for the P(+)-VX stereoisomer. Although these numbers cannot be construed as true LD₅₀ values, they are comparable to the LD₅₀ values identified by Wright et al. (2017).

Table 3. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for LD₅₀ Studies

Guinea	Exposure	Dose	Sample	Lavellous Lapos	· /	V	K-G																			
Pig No.	Date	(mg/kg)	Time	(ng/g)																						
1191100	2000	(8/8/	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain																	
11]		16	367.6197	519.5762	820.5367	3,512.6174	4,128.8413	350.9682																	
12			15	215.1405	461.4225	674.4045	2,257.0621	5,888.4256	376.4869																	
13			15	447.3438	510.7501	626.7815	1,704.0373	3,107.1309	571.0728																	
14			13	299.8104	542.5614	746.6021	1,782.4065	2,994.5421	578.2536																	
16	1	1.0*	17	292.8682	404.1579	599.5910	2,306.9841	1,754.6218	522.8988																	
17		1.0	13	388.9109	696.2036	862.6254	6,365.2048	1,156.4854	802.8230																	
19			18	269.63	357.0062	305.9527	3,382.8852	3,162.8028	668.2355																	
20			13	288.6617	578.7843	738.5032	1,526.3359	1,663.6786	488.7705																	
21			16	216.6482	616.7007	1,008.0593	775.3744	1,325.6986	643.1431																	
22			16	440.6334	483.3809	782.4782	2,358.8727	2,943.0750	320.3121																	
15] [3.8512	4.5828	17.7468	154.5606	134.2647	2.7563																	
23				3.5665	2.5943	10.2699	107.5471	449.6310	4.2018																	
24	27 Jan 2016			3.7626	2.1769	9.4407	157.8269	317.7781	3.8584																	
25				4.0261	2.2869	12.8848	139.9483	90.0866	3.9150																	
26		0.14	0.14	0.14	0.1*	1 440	3.0958	3.8437	23.4856	421.1333	699.7246	5.0165														
27		0.1	1,440	3.4843	4.0521	37.6090	103.8152	609.2639	4.2143																	
28			ı														1				3.1834	2.9490	17.2232	242.5674	540.1696	3.4539
29]			3.6057	3.0589	13.8423	154.4708	606.1807	3.1835																	
30				3.5472	3.4553	11.9384	250.6466	441.2876	4.7191																	
31				3.3949	2.8644	11.8563	130.3585	344.4878	3.8872																	
33] [30	135.7689	213.0770	349.4411	1,684.4153	2,226.6943	138.0405																	
34]		31	123.5151	154.7620	277.4501	4,308.0255	8,963.2723	111.8013																	
35		0.5*	27	105.6681	158.0451	320.8774	2,865.1498	4,345.6519	120.4120																	
37]		35	135.7413	137.1609	250.6580	2,477.7080	2,561.2895	93.5203																	
38]		33	103.6892	112.8116	232.7419	2,702.8533	2,958.5497	94.8317																	

^{*2-}Propanol was used as the solvent for these dosings.

Table 3. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for LD50 Studies (Continued)

Guinea	Exposure	Dose	Sample Time	VX-G (ng/g)									
Pig No.	Pig No. Date (m	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain				
39			30	98.2962	132.6429	303.0177	1,888.2508	4,931.0791	107.3883				
40			32	95.5986	130.0244	265.5998	3,413.2856	2,473.6363	88.2459				
41	27 Jan 2016	0.5*	0.5*	0.5*	33	128.5214	132.5289	228.2742	2,662.1422	1,674.2905	106.7796		
43			28	No sample	154.6299	259.3492	2,718.8621	8,030.0944	125.5538				
44			29	117.1229	132.4259	129.8050	2,075.2998	4,815.7786	114.4125				
45			40	90.4188	100.1918	170.6575	2,448.9813	1,911.9372	69.4550				
46		0.3*	0.3*	48	74.1967	61.8709	114.7646	1,263.0950	1,525.4601	37.5772			
48				88	37.6565	35.8046	105.0481	2,951.7492	5,814.3531	18.9562			
49							46	75.4668	89.1963	184.2249	2,645.4985	1,452.8385	82.0056
50				45	71.9854	75.8710	181.1715	2,374.3413	1,485.9552	56.7199			
51			42	68.5553	85.9407	179.6666	1,676.2200	2,549.1703	56.7042				
52			45	67.6037	64.7845	161.8756	1,622.5549	3,533.9418	39.2166				
53			39	98.1148	85.9370	184.3899	1,681.5725	960.4025	52.9353				
55	28 Jan 2016		40	89.2016	111.2657	227.5087	1,183.8118	2,088.2020	86.1643				
57				4.1804	10.0335	72.9183	908.6338	1,132.4542	6.0293				
58				3.7631	12.2198	38.9625	944.2015	1,070.4922	5.0376				
59				5.5608	13.9915	55.3750	512.5140	285.4990	5.5768				
60		0.15*	1,440	4.6636	12.0834	61.6722	938.4950	336.5950	6.9748				
61		0.13	1,440	7.7892	10.4946	48.4688	521.0721	3,228.3681	6.1962				
62				6.1386	7.5356	38.0267	366.3186	575.4307	5.7104				
64				4.5203	7.6356	18.2784	146.0636	56.2294	3.8817				
36				4.3186	5.7927	24.8023	441.0047	1,444.0826	4.6796				
1			53	94.8408	75.4575	141.1850	2,096.9689	4,677.3025	72.9816				
3	3 Feb 2016	0.21*	73	63.2023	49.5303	147.1370	1,712.1490	8,233.5432	26.5190				
4	3 Feb 2010	0.21	76	69.1603	43.7102	184.1990	2,652.6278	7,584.2297	27.2607				
5			78	75.0011	45.8545	108.3745	2,251.4199	1,479.2198	21.9598				

^{*2-}Propanol was used as the solvent for these dosings.

Table 3. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for LD₅₀ Studies (Continued)

Guinea	Exposure	Dose	Sample Time				K-G g/g)				
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
6			44	109.4726	80.4968	174.9116	1,206.3224	1,677.4982	60.8008		
7			51	90.9247	80.3956	123.3876	1,228.2390	730.0775	53.2763		
47	3 Feb 2016	0.21*	55	71.7622	66.5633	384.8760	2,767.4731	3,913.3086	41.8769		
54			54	58.9190	57.5608	128.6920	1,780.3813	1,687.4593	41.7599		
65	7		55	51.0443	90.5181	205.8334	2,900.4081	4,889.6125	70.7391		
221				4.8167	57.4764	154.6742	2,445.1692	1,355.4395	17.0329		
223				4.0630	94.2406	366.7278	2,268.8692	2,913.9389	23.2031		
224				22.3542	226.1358	796.0360	1,600.4834	2,710.0788	56.0006		
225				4.7670	221258	91.9891	1,116.0255	979.9203	9.2807		
226	5 1 2016	0.175	1,440	4.5206	26.2746	89.8677	802.5146	448.4934	11.3057		
227	5 Apr 2016	0.175		4.7093	18.4959	91.3592	1,201.2346	148.1926	9.6939		
228				2.8895	31.7032	87.5128	849.1210	569.1677	5.6206		
230				5.7890	1.8069	6.5525	301.4764	684.5524	7.1732		
231				14.7563	30.2467	231.1992	572.7788	985.0512	14.2989		
233				3.4676	11.5993	56.2404	80.9442	165.1171	10.7063		
234				2.3783	1.0330	5.5495	6.1166	14.5521	1.7457		
235						2.4416	0.9642	9.8437	2.6155	39.4614	1.2475
236				3.2765	1.2011	3.8910	1.8321	20.6584	1.1668		
237				2.3516	1.9806	4.4907	1.7497	14.1896	1.1103		
238	27 Apr 2016	0.062	1,440	2.0033	1.0363	4.8924	3.0372	64.7529	1.4033		
239				2.2244	1.0476	4.8973	4.5923	51.8611	1.0602		
240				2.8251	0.9896	3.7095	2.0536	68.3028	1.2595		
241				2.8519	1.0917	3.0039	2.6498	10.3481	0.8283		
242				2.6258	0.8795	3.5483	4.3647	31.0824	0.9999		
243	28 Apr 2016	0.28	35	108.5499	103.8073	450.0490	1,128.7090	5,588.0942	105.8788		
244	20 Apr 2010	0.20	40	77.3273	106.1143	175.8353	794.4640	3,959.6067	60.7724		

^{*2-}Propanol was used as the solvent for these dosings.

Table 3. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for LD_{50} Studies (Continued)

Guinea	Exposure	Dose	Sample Time	VX-G (ng/g)						
Pig No.	Date	Date (mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
245			86	48.5002	107.2460	220.2025	597.8624	2,850.0274	26.4602	
246	1		34	124.0270	101.2686	182.1559	1,470.4289	2,156.0354	94.4460	
247]		37	93.1963	85.5351	184.9963	4,272.7685	5,109.9987	58.3185	
248	28 Apr 2016	0.28	52	57.2503	51.5312	124.3433	795.9576	915.9505	50.3942	
249]		96	40.0441	19.0486	78.9444	1,613.5982	1,450.5470	14.4586	
250]		103	28.0954	17.4541	121.5443	1,014.3769	2,279.3179	18.1765	
251			98	28.4870	16.7098	69.0781	1,332.9761	1,500.6994	18.0519	
273				5.9890	13.6205	30.9093	91.6437	100.7163	7.1782	
274	21 Jun 2016	0.21	1,440	4.5176	2.1578	12.2321	71.6553	139.9559	7.9459	
275				8.2419	20.1642	109.5028	86.2091	368.9131	10.1986	
276			1,440	4.2281	15.8002	40.2624	104.3492	48.9155	8.2499	
277				3.9404	4.1862	28.7462	30.3145	54.7813	4.8221	
278		0.21		3.0472	35.7713	57.0586	311.5609	618.7852	11.0895	
279				5.3044	24.3384	80.2497	173.3943	478.4172	10.2999	
280				4.8842	7.0342	23.5295	77.8334	161.4437	7.1219	
281				3.6383	6.8727	15.4113	372.6103	115.5208	6.0203	
283				3.7483	19.9105	43.7904	442.0844	399.0313	11.4181	
284	21 Jun 2016			4.5892	11.0285	43.0022	57.3478	233.1606	8.2986	
285	21 Juli 2010		1,440	4.6196	15.5057	39.5968	148.7849	511.8685	9.6933	
287			1,440	3.7718	19.1143	65.6033	165.4153	69.9246	9.8110	
288]	0.245		4.0112	10.6957	24.8123	808.7154	567.6715	8.8511	
289		0.243		3.8825	8.5966	31.6284	135.9029	192.6893	6.7034	
290			<1,440			No sa	mples			
291				3.6427	14.2433	49.2973	177.2697	61.9793	9.6805	
292			1,440	3.3932	26.2629	61.2158	54.7325	67.6390	6.9643	
293				4.1152	15.4529	46.4755	245.6619	156.3586	8.3449	

<1440, Animal died prior to 24 h.

11

Table 3. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for LD₅₀ Studies (Continued)

Guinea Exposure Pig No. Date	_	Dose	Sample Time	VX-G (ng/g)						
	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
294			1,440	6.2277	17.0478	47.2814	46.3610	66.0524	6.5113	
297		0.26	1,440	4.7024	8.8271	9.8100	49.2839	118.2668	6.2357	
298			1,440	3.5244	7.4176	15.3535	46.9616	83.2553	5.1433	
299	22 Jun 2016		1,440	2.3183	10.9795	19.0715	61.7500	70.5455	6.1543	
300			1,440	2.8816	8.9481	20.3020	41.1902	28.6108	5.2837	
301			1,440	5.8271	9.2858	18.7115	43.4706	62.0381	8.8871	
302			1,398	68.7272	35.0618	18.7101	44.6110	88.3184	4.7413	

Table 4. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for LD₅₀ Studies

Guinea	Exposure	Dose	Sample Time	VX-G (ng/g)						
Pig No. Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
68			30	5.3347	3.5900	6.8047	9.1295	22.9099	6.2092	
79			29	5.0861	4.3729	8.8116	8.5498	27.0105	9.9515	
85			37	4.0854	2.6730	4.8506	6.0027	8.9081	8.7434	
92			19	5.7929	5.1043	7.7305	3.4509	16.2395	10.0177	
103	1 Mar 2016	Mar 2016 3.5 ^a	29	4.8005	3.9759	9.3099	8.3823	13.2611	6.6656	
108			45	5.4250	2.8958	5.6320	8.3199	9.3960	8.9735	
114			23	4.9100	3.6724	8.3468	6.3332	10.2269	9.0837	
120			24	4.8573	4.3357	7.6110	4.8612	20.8456	3.9985	
125			37	4.7372	3.3865	5.9785	9.4440	14.1312	9.6679	
66	7 Mar 2016 2.2 ^{a,b}	1 440	0.7480	0.3604	0.6018	0.3432	0.9463	2.6316		
67	7 Mar 2016	2.2.,.	1,440	1.0077	0.6987	1.3422	0.4304	1.1254	1.5899	

^a2-Propanol was used as the solvent for these dosings.

^bFollowing animal exposures, NMR analyses of neat agent used to prepare these dosings indicated only 81–83% purity.

Table 4. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for LD₅₀ Studies (Continued)

Guinea	Exposure	Dose	Sample Time	LAPOSUIC to I		VX	(-G g/g)		
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain
69				0.5695	0.3114	0.4977	0.2877	BDL	2.5807
70				0.6527	0.3346	0.4340	0.2911	0.5929	1.0903
71				0.6752	0.4133	0.5562	0.3694	BDL	2.5455
72		$2.2^{a,b}$	1 440	0.7433	0.3135	0.3577	0.2930	BDL	2.5874
74		2.2	1,440	0.7097	0.3419	0.6918	0.2664	BDL	2.0434
76				0.6999	0.2941	0.4117	0.1426	1.4507	1.8837
81				0.8296	0.3765	0.6580	0.3941	BDL	3.3260
82	7 Mar 2016			0.7366	0.3782	0.7192	0.2081	BDL	2.9569
84	7 Mar 2016		1,440	0.9351	0.3934	1.0883	0.4157	0.7245	3.5956
86			1,440	0.7629	0.4074	0.7293	0.5247	BDL	3.0410
90		3.0 ^{a,b}	67	3.5862	1.4655	2.8489	2.7338	14.4339	5.992
91			1,440	0.7135	0.5821	0.8580	1.4001	1.0534	6.1766
94			1,440	1.0824	0.8750	1.2565	0.8584	1.2254	3.6520
131				0.8807	0.5714	0.6971	0.4373	2.3775	3.3456
134		$3.0^{a,b}$	1,440	0.9230	0.3735	0.6085	0.4757	BDL	3.1044
140				0.7777	0.4407	0.8644	0.4868	BDL	4.2573
77				0.6272	0.3436	0.6666	0.4165	BDL	3.4023
88				1.1030	0.2793	0.3418	0.5630	BDL	3.1018
93				0.8436	0.2931	0.3102	0.5717	BDL	2.1655
100				0.8473	0.4823	0.6859	0.9196	2.3387	4.7672
101	22 Mar 2016	3.0	1,440	1.0259	0.4300	0.7899	0.6471	0.9119	4.6040
102	22 Mar 2016			0.9023	0.4422	0.5599	0.4348	1.3578	3.0996
128				0.9120	0.2890	0.7004	0.4942	0.6432	3.7619
129				0.8065	0.3198	0.3425	0.6894	0.9009	2.9097
130				0.7579	0.4674	0.2777	0.5255	0.8636	3.2750
95		5.0	35	4.6482	2.8559	6.0981	10.5137	28.6588	8.0786

^a2-Propanol was used as the solvent for these dosings.

^bFollowing animal exposures, NMR analyses of neat agent used to prepare these dosings indicated only 81–83% purity. BDL, below detection limit (<0.05 ng/g).

Table 4. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for LD₅₀ Studies (Continued)

Guinea	Exposure	Dose	Sample Time				(-G g/g)									
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain							
96			181	2.9662	1.4668	3.5134	6.4265	9.4069	8.7913							
105			113	3.4499	1.5913	4.9658	7.9365	12.2631	5.4513							
109			134	2.9156	1.5205	2.4912	10.5164	18.0524	6.2803							
113		5.0	1,440	0.9616	0.4036	0.7960	1.6326	5.3464	6.6717							
132		3.0	61	4.1542	2.1711	4.4736	3.7455	13.8886	9.7702							
133			172	3.3345	1.8309	4.2097	10.4734	12.8433	11.0950							
137			85	2.7058	1.7109	4.2721	5.1679	35.0365	10.9229							
138			36	3.7383	2.7063	4.9044	5.7604	19.0318	13.2501							
97	22 Mar 2016		1,440	1.0939	0.2897	0.4408	0.3764	0.4261	4.0612							
115	-22 Mar 2016		1,440	0.6201	0.3895	0.4856	0.7341	1.0280	5.9691							
116		4.0	146	2.7199	1.0828	2.4762	3.4731	10.5987	7.2884							
139			1,440	4.6341	1.0809	1.4369	1.9357	27.8021	6.9294							
141			146	2.6236	1.2311	4.8061	3.3888	5.7368	3.8039							
147			1,440	2.9745	2.4337	4.0501	8.8499	51.7553	11.9443							
148			51	0.9034	0.9463	0.7575	0.8917	7.1966	0.1800							
211				0.7144	1.6111	3.2241	3.1373	61.8846	6.6070							
212		4.2	4.2	4.2	4.2	4.2	4.2	4.2	4.2	1,440	0.9331	0.3575	0.4538	0.7813	16.1264	4.5767
213				0.9604	0.3067	1.8147	1.8192	23.5208	3.7092							
214			1,440	0.6444	0.4597	4.3666	0.8603	8.6556	2.5417							
216	22 Mar 2016 4.2	242	2.5964	1.2692	4.3196	7.7641	22.7743	7.0768								
217		4.2	1,440	0.9852	0.4966	2.5123	1.8837	3.0178	4.0519							
218		4.2	1,440	0.7781	0.4648	2.0600	2.6934	2.3115	5.4677							
219			1,440	1.2112	1.5973	3.4611	0.5725	0.7330	1.6194							
220			1,440	0.8676	0.9498	3.9424	1.2286	4.1802	5.2600							

Table 5. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for LD₅₀ Studies

Guinea	Exposure	Dose	Sample				(-G		
Pig No.	Date	(μg/kg)	Time			(ng		T	
	2	(F-8/8/	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain
143			17	15.3924	16.8396	27.8820	23.0935	116.4108	28.8881
146			15	11.0238	14.5599	40.4316	25.7339	277.5318	22.1242
149			16	10.3654	13.9037	27.3427	33.4258	43.4299	22.3774
150	7 Mar 2016	7.0*	13	19.6014	19.8311	58.2037	17.5263	68.7084	19.9897
152	/ Wiai 2010	7.0	12	12.9509	18.3767	31.7553	23.6811	251.2938	35.1998
166			19	11.2233	10.2837	22.3453	85.5740	268.0517	15.9662
167			13	16.0250	19.4756	44.4074	12.0133	243.8794	24.5218
168			18	10.3003	14.7997	37.0817	64.2143	433.8273	22.2755
78			42	5.5961	4.5814	10.5038	31.7573	103.7902	11.0781
83			33	7.0021	5.2740	16.4919	31.9507	95.6407	10.0623
89		7.0	47	7.1829	5.3989	11.4316	58.2181	227.8316	14.1338
107			1,440	6.4610	5.5808	11.5955	43.9200	164.3720	12.0488
123			1,440	7.2369	4.6701	9.8530	35.6108	269.1729	12.6398
111	1			1.2676	0.5826	1.9367	7.6974	16.3417	3.5263
112			1,440	1.0183	0.5411	0.7833	3.5242	2.2788	2.8692
136	10 Man 2016			1.2161	0.7971	1.8774	5.4279	17.7490	2.3890
169	10 Mar 2016	5.0	<1,440	1.4283	0.6255	1.1684	6.9638	18.6265	4.2002
170				0.9755	0.6583	3.0474	12.6399	13.8405	7.5253
171			1,440	1.1947	0.5716	1.6880	11.7517	27.3810	4.2695
172				0.4199	0.6818	3.4606	6.5627	5.5444	4.2980
75] [1,440	1.2171	0.8112	2.3843	8.8715	23.3525	5.9716
119]	6.0	28	3.7418	3.4580	5.8318	40.6825	143.6869	10.5019
176	6.	6.0	1,440	1.3215	0.5127	1.8625	10.0747	19.1393	3.8866
178			1,440	1.3517	0.6316	1.6762	12.9665	22.5888	4.5305
153	23 Mar 2016	4.8	171	6.3430	0.2586	16.5318	40.3196	42.7763	9.0433

^{*2-}Propanol was used as the solvent for this dosing. <1440, Animal died prior to 24 h.

Table 5. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for LD_{50} Studies (Continued)

Guinea	Exposure	Dose	Sample Time		VX-G (ng/g)						
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
154		4.8	1,440	2.5391	1.1059	2.6704	8.5782	11.4143	5.0702		
155] [1.8622	1.4041	2.8306	7.3260	8.8541	2.3452		
156]	4.1	1 440	2.0305	1.4473	3.6337	4.3505	7.7146	2.3278		
157			1,440	2.0672	0.7752	2.0898	6.8690	24.8946	4.6872		
158				1.7007	0.9460	3.0389	5.7417	21.5219	4.2866		
127	23 Mar 2016			1.9959	1.3496	2.5106	4.3980	4.3915	3.7892		
161	23 Mar 2010			2.0462	1.6088	7.9095	15.9133	23.9412	4.3646		
163				1.7888	2.0293	1.7410	4.0250	8.4988	2.5042		
164	3.7	3.7	1,440	2.0050	1.1825	2.8618	8.6070	21.6200	3.7304		
173				1.8237	1.0409	2.9634	8.8657	14.6460	5.2388		
175				1.6983	1.4677	4.7437	5.8650	5.2279	4.8160		
179				1.7956	1.3872	2.3457	6.8139	33.0112	3.2593		




Figure 2. VX-G concentration in blood as a function of the dose of individual VX stereoisomers or the racemic mixture.

3.3 Pharmacokinetic Studies

Table 6 summarizes the data from the VX-G assays of whole blood, tissues, and organs that were obtained from guinea pigs after intravenous exposure to various doses of P(+)-VX and collection of serial blood samples. Tissues and organs were only harvested after death or at the time of euthanasia. Tables 7 and 8 show similar results from intravenous exposures to various doses of P(-)-VX and a racemic mixture of VX, respectively.

Figure 3 is a plot of the VX-G concentration in blood as a function of time for three different doses of P(+)-VX. Each point represents an average (n = 4 animals) of the individual values that were determined at the time of sampling for each dose prepared in saline. Error bars show the standard deviation of each average determination. Curves were fitted using a two-parameter pharmacokinetic model of the form $y = a \times e^{-b(x)}$, where x and y represent the corresponding axis values on Figure 3, a is the VX-G concentration at time 0, and b is the elimination rate constant (which is the same as k in Table 9). This assumes that after administration of an intravenous bolus of VX, the transfer of VX from the blood (as measured by the VX-G concentration) follows first-order kinetics. Similar data analyses were performed from serial samplings after intravenous exposures to various doses of P(-)-VX and a racemic mixture

of VX. Table 9 summarizes the elimination rate constants (k, in inverse minutes) for the various doses of VX stereoisomers and the racemic mixture as determined by the slope of the pharmacokinetic models. The half-life ($t_{1/2}$) of VX (as VX-G) remaining can be calculated as $t_{1/2} = \ln(2)/k$, as is also shown. The data revealed that a much slower clearance rate (15–20 times slower) occurred for the P(–)-VX stereoisomer and the racemic mixture, which correlates with the higher toxicity of the P(–)-VX stereoisomer and the racemic mixture.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies

Guinea	Exposure	Dose	Sample Time		71 (+)- V X 101 1 1	VX (ng	C-G					
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
205			0			No sample	es received					
			13	47.3661	90.4690	161.9064	1,479.0200	2,744.6300	82.4140			
			0									
			1	ISx								
			10	153								
	12 Apr 2016	0.175	30									
207			40] ,		No	o samples receive	ed				
207			50	Sample								
			60	assays were lost during								
			180	preparation								
			360	• •	_							
			1,440	1.8742	8.8404	16.0046	553.1206	302.1590	7.3588			
			0	BDL								
			1	52.8231								
			10	32.6265								
			20	27.6206								
			30	20.3389		No	o samples receive	ed				
252	1 Jun 2016	0.175	40	17.9902		140	samples receive	cu				
			50	15.2675								
			60	12.9064								
			180	5.7845								
			360	5.0676			<u> </u>	,				
			1,440	3.0659	15.7813	36.9404	73.8697	66.5627	6.5117			

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time	posure to F(+)- v		VX	Κ-G g/g)		
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain
			0	BDL				1	
			1	ISx					
			10	33.2948					
254		0.175	20	33.1621		N	o samples receiv	ved	
234		0.173	30	ISx					
			40	31.8474					
			50	33.0211					
			64	16.7881	23.3540	40.7594	2,726.6284	1,107.5617	23.2660
			0	BDL					
			1	ISx					
			10	157					
			20	46.3713		N	No samples received		
268			30	32.5281		11	o samples receiv	red	
	1 Jun 2016		40						
			50	ISx					
			60						<u>.</u>
	<u></u>		203	20.6665	7.2703	53.0362	776.0446	826.5385	9.0034
		0.28	0	BDL					
			1	107.9733					
			10	52.7423					
			20	40.5516					
270			30	39.0246		N	o samples receiv	ved	
2.0			40	28.4175		1,	o samples receiv		
			50	25.0706					
			60	19.1688					
			180	14.5818					
			360	11.6060					

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose (mg/kg)	Sample Time				X-G g/g)			
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
			0	BDL						
			1	103.8394		N				
			10	59.9017		IN	o samples receiv	ed		
271			20	58.2990						
			30			N	o commiss massiv	a d		
			40	ISx		IN	o samples receiv	ed		
			63	1	30.0360	93.7820	516.6617	212.8916	28.6461	
]		0	BDL						
	1 Jun 2016	0.28	1	75.9267						
			10	34.0414						
			20	25.3866						
272			30	22.6829		N				
272			40	19.7902		IN	o samples receiv	ed		
			50	14.8970						
			60	13.2719						
			180	6.1499						
			360	4.6188						
			1	33.4396						
			10	27.1744						
262	12 Jun 2016	0.21	20	22.8026		N	o samples receiv	red		
262	13 Jun 2016	0.21	50	28.1692	8.1692					
			60	IC	IC.,					
			153	ISx	X 31.8201 57.3465 516.5460 1,425.8684 11.650					

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time	posure to F (+)-V		VX	K-G g/g)			
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
			1	63.9592		•	•			
			10	34.1746						
			20	26.5895						
			30	21.8265						
266			40	19.5429		No	Samples Receiv	ved		
200			50	16.6621						
			60	16.4341						
	13 Jun 2016	0.21	180	11.2967						
	13 Juli 2010	0.21	360	Samples lost	2.9776 25.1167 323.5490 518.2482 6.661					
			1,440	Samples lost						
269			10	ISx		N	o samples receiv	ad		
209			30	29.3581		11	o samples receiv	eu		
			40	ISx						
269			50	19.2091		N	o samples receiv	ed		
209			60	24.4821						
			239	BDL	7.4591	39.3587	248.0948	467.3637	8.8173	
			0	BDL						
			10	68.4203						
			20	ISx						
			30	42.1468	No samples received					
306	29 Jun 2016	0.26	40	ISx		14	o samples receiv	Cu		
			50	31.8987						
			60	25.8275	5.8275 ISx					
			180							
			335	20.6178	7.8376	37.7228	408.5275	229.3400	6.7097	

BDL, below detection limit (<0.05 ng/g).
ISx, an insufficient amount of sample was received for assay.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample				(-G g/g)					
Pig No.	Date	(mg/kg)	Time (min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	10.1510								
			1	11.4679								
			10	22.2821								
			20	35.5853								
			30	32.9902		NT.	No samples received					
308			40	25.1858		INC	o samples receiv	ea				
			50	24.6068								
			60	47.2548								
		180	70.9735	1								
		360	0.6111	1								
	29 Jun 2016	0.26	1,440	4.3191	1.8057	14.6519	106.9136	144.5937	4.338			
			0	BDL								
			1	118.6662								
			10	50.5825								
			20	38.3982		No	o samples receiv	ed				
200			30	29.6539								
309	309		40	25.0645								
			50	22.9471								
			60	19.4047								
			180	7.8441	No samples received							
			360	7.6641								

BDL, below detection limit (<0.05 ng/g).

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time				K-G g/g)				
No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
			0	BDL							
			1	111.6205							
			10	49.2589							
			20	35.8388							
310		0.26	30	30.1986		No	o samples receiv	ed			
310		0.20	40	21.9792							
			50	19.8471							
			60	16.0746							
			360	7.1977							
	20 I 2017		1,440	9.2727	36.9375	233.8527	121.2871	73.4495	6.7068		
	29 Jun 2017		0	BDL							
			1	ISx							
			10	31.5900							
			20	25.6230							
			30	23.5765		N	o samples receiv	ad			
311		0.245	40	18.4448		111	o samples receiv	eu			
			50	17.3423	7.3423						
			60	13.7759							
			180	5.6444							
			360	5.4667							
			1,440	3.8136	1.6437	5.3239	65.1485	37.7119	2.6862		

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time	posure to F(+)-v		VX (ng	[-G			
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
			1	119.8871						
			10	51.6296						
			20	43.0279						
303			30	34.3514		No	samples receiv	ed		
			40	28.5456						
			50	27.5613						
			60	ISx						
			180	9.3291		N	a complac racaiy	rad		
303			360	6.5514	No samples received 1 7756					
			1,440	4.2221	1.7756	12.7932	78.9552	51.0077	5.6433	
			0	BDL						
			1	110.2506						
	6 Jul 2016	0.245	10	53.0537						
	0 341 2010	0.243	20	40.8154						
304			30	35.4353		N	samples receiv	rad		
304			40	30.8270		140	o samples receiv	eu		
			50	26.1509						
			60	ISx						
			180	14.4464						
			360	7.8400						
			1	85.0553						
			10	45.3622						
305			20	40.0229		NI	samples receiv	ed		
303			30	33.2793		110	samples lecely	eu		
			40	32.4854						
			50	34.9095						

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

25

Table 6. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(+)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose (mg/kg)	Sample Time	000010 10 1 (1)			X-G g/g)			
Pig No.	Date	(mg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
			60	ISx						
305		0.245	180	8.5876		N	lo samples receiv	ed		
303		0.243	360	4.9088						
			1,440	3.7082	2.0069	8.8287	74.9796	53.2529	5.9840	
			1	99.6083						
			10	53.9137						
			20	38.4248						
307		0.21	30	37.0028		N				
307		0.21	40	26.5192		IN	lo samples receiv	ed		
			50	26.2970	26.2970					
			60	18.9385						
	6 Jul 2016		180	8.9507						
307		0.21	1,440	7.2095		N	lo samples receiv	ed		
			1							
			10							
			20							
			30							
212		0.175	40	ISx	ISx No samples received					
312		0.175	50							
			60							
			180							
			360							
			1,440	4.7165	8.8235	18.0489	40.4379	33.0470	6.3909	

ISx, an insufficient amount of sample was received for assay.

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies

Guinea Pig No.	Exposure Date	Dose (µg/kg)	Sample Time (min)	VX-G (ng/g)							
				Whole Blood	Heart	Lung	Liver	Kidney	Brain		
	5 Apr 2016	4.2	0	BDL							
			1	ISx							
			10	2.3125							
			20	3.4517							
			30	2.6101							
188			40	2.5747	No samples received						
			50	2.8179							
			60	2.6433							
			180	2.2861							
			360	1.9649							
			1,440	0.7280	0.3518	0.4533	0.9912	2.8583	3.8347		
			0	BDL							
			1	2.3679							
			10	2.2212							
			20	2.3780							
189			30	2.4598	No samples received						
109			40	2.7345	ino samples received						
			50	2.6416							
			60	2.5074							
			180	2.2846							
			360	1.7257							

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea Pig No.	Exposure Date	Dose (µg/kg)	Sample Time (min)	VX-G (ng/g)							
				Whole Blood	Heart	Lung	Liver	Kidney	Brain		
	5 Apr 2016	4.2	0	111%*							
			1	2.6712							
			10	3.1222							
			20	3.3248							
			30	3.1491	No samples received						
192			40	2.2934							
			50	2.8705							
			60	3.0374							
			180	2.5376							
			360	2.1062							
			1,440	0.8312	0.3891	0.4352	0.5124	1.3448	3.0172		
194			0	ISx							
			20	ISx	No samples received						
			40	ISx							
			50	3.1004							
			60	2.9951							
			180	2.7689							
			360	2.1006							
			1,440	0.8987	0.4244	0.5582	0.9501	0.8233	6.3418		

^{*}This sample was spiked with VX-G to determine percentage recovered from assay.

(continued)

ISx, an insufficient amount of sample was received for assay.

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea Pig No.	Exposure Date	Dose (µg/kg)	Sample Time (min)	VX-G (ng/g)							
				Whole Blood	Heart	Lung	Liver	Kidney	Brain		
209	5 Apr 2016	4.0	0	BDL							
			1	1.1506							
			10	2.4687	No samples received						
			20	2.8329							
			30	2.8066							
			40	2.9751							
			50	2.6572							
			60	2.7573							
			180	2.4735							
			360	2.0305							
			1,440	0.7594	0.3252	0.3970	0.7829	1.1482	4.2660		
	31 May 2016	3.0	0	BDL	No samples received						
253			1	2.1660							
			10	2.0163							
			20	2.0864							
			30	1.9685							
			40	1.6704	No samples received						
			50	1.8538							
			60	1.9363							
			180	1.5611							
			360	1.4948							
			1,440	0.3707	0.2453	0.6178	0.4852	1.5160	2.1775		

BDL, below detection limit (<0.05 ng/g).

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time				X-G g/g)					
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL			1					
			20	2.7002								
			30	2.2413			No samples received					
			40	2.4665		N						
255			50	1.8995		18						
			60	1.7895								
			180	1.8420								
			360	1.5622			0.3831 0.2120 0.8171 2.663					
			1,440	0.5641	0.2625	0.3831						
		3.0	0	BDL								
		3.0	1	3.4256								
			10	2.6099								
			20	2.7143								
	31 May 2016		30	2.8300		N	o samples receiv	ad				
256			40	2.7376		11	o samples feceiv	cu				
			50	2.4541								
			60	2.6456								
			180	2.2713								
			360	1.8868								
			1,440	0.5318	0.2634	0.4126	0.2845	8.1476	3.2098			
			0	BDL								
			1	3.0384								
			10	3.0757								
257		4.0	20	3.6126	No samples received							
			30	3.2951			_					
			40	3.3537								
			50	3.2887								

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose (µg/kg)	Sample Time		, 11 101 1 11	V	K-G g/g)				
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
			60	3.4274							
257			180	2.7760		N	o samples receiv	red			
	_		360	2.0940							
			0	BDL							
			1	2.9114							
			10	2.6274	No samples received						
			20	3.0077							
258			30	3.0697							
230			40	3.2345							
			50	2.8910							
			60	2.8104							
	31 May 2016	4.0	180	ISx							
	31 Way 2010	4.0	360	1.4351							
			0	BDL							
			1	2.9977							
			10	2.8593							
			20	2.9727							
			30	3.0966	_	N	o samples receiv	red			
259			40	2.8144	<u> </u>						
			50	2.8352							
			60	2.6607	_						
			180	2.4734							
			360	1.9287							
			1,440	0.6568	0.3279						

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample			V	X-G g/g)				
Pig No.	Date	(µg/kg)	Time (min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
			0	BDL							
			1	3.1315							
			10	3.1833							
			20	3.5888							
			30	3.6730		No samples received					
260			40	1.4309							
			50	3.7682							
			60	3.8095	_						
			180	2.9990							
			360	2.5796							
	31 May 2016	5.0	1,440	0.7426	0.4592	0.6490	0.6488	0.5637	3.2290		
			0	BDL	_						
			1	3.1648							
			10	3.7401	_						
			20	3.9175							
261			30	3.5355		N	o samples receiv	ed			
201			40	4.0766		11	o samples receiv	cu			
			50	3.7690							
			60	3.6321							
			180	2.6953							
			360	1.9450							

BDL, below detection limit (<0.05 ng/g).

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time				X-G g/g)						
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain				
			0	BDL									
			1	3.5571									
			10	3.0994									
			20	3.0729									
264		5.0	30	3.3483		N	o samples receiv	ad					
204		5.0	40	3.7155	55								
			50	3.4073									
			60	3.2426									
			180	2.5005									
	_		360	2.2232									
	31 May 2016		0	BDL									
			1	2.0273									
			10	2.1867									
			20	2.3503									
265			30	2.3126		N	o samples receive	ed					
		3.0	40	2.6012									
			50	2.3676									
			60	2.3346									
	_		180	2.0540									
265			360	1.7699	ļ		o samples receiv						
203			1,440	0.5594	0.4004	0.3819	0.1220	4.6793	3.7362				

Table 7. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to P(–)-VX for Pharmacokinetic Studies (Continued)

Guinea	Exposure	Dose	Sample Time				K-G g/g)					
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL								
			1	3.2516								
			10	3.3887								
			20	3.4344		N	o samples receiv	ad				
267	31 May 2016	5.0	30	3.2779		IN	o samples receiv	eu				
			40	3.3340								
			50	3.3395	395							
			60	2.9954								
			167	1.8184	1.3027	2.5009	2.6029	24.6279	5.9398			

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies

Guinea Big No	Exposure Date	Dose (ug/lrg)	Sample Time				K-G g/g)					
Pig No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL								
			1	5.6508								
			10	4.9878								
			20	4.0480	No samples received							
			30	4.9290								
180			40	4.2287								
			50	4.1203								
			60	4.2309								
			180	3.8336								
			360	3.1521								
	23 Mar 2016	4.8	1,440	1.6439	0.8790	1.6369	12.5498	26.5548	3.7624			
	23 Wai 2010	4.0	0	BDL								
			1	5.4642								
			10	5.8777								
			20	5.1333								
			30	5.7275		N	o samples receiv	red				
182			40	5.3119		14	o samples receiv	cu				
			50	4.7095								
			60	5.2397								
			180	4.6352	63							
			360	3.5663								
			1,440	1.6162								

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	dire to Raceinic		VX	K-G g/g)	,				
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL			1	•	•			
			1	5.8179								
			10	5.1014								
			20	5.5445								
			30	5.5172		No samples received						
183			40	5.3774	No samples received							
			50	5.1176								
			60	5.1797								
			180	4.4672								
			360	3.7489								
	23 Mar 2016	4.8	1,440	1.8634	1.1789	1.8316	11.0851	26.2179	6.3050			
	23 Wai 2010	4.0	0	BDL								
			1	5.9456								
			10	5.2385								
			20	4.8052								
			30	4.4030		N	o samples receiv	red.				
185			40	4.4884		11	o samples receiv	eu				
			50	4.2926								
			60	3.9518								
			180	3.5226								
	360 2.9443											
	1: 1: 1: 1:		1,440	1.4662	.4662 0.6950 1.4151 15.0011 23.3111 4.7752							

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	The real real real real real real real rea		VX	K-G g/g)	,				
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL								
			1	3.9543								
			10	3.4025								
			20									
			30		No samples received							
181			40	ISx								
			50	157								
			60									
			180									
			360	2.2645								
	30 Mar 2016	7.0	1,440	1.0349	0.9668	2.1490	9.4032	171.7250	3.5190			
	50 1/1 41 2010	7.0	0									
			1	ISx								
			10									
			20	4.2198								
			30			N	o samples receiv	red				
184			40		-							
			50	ISx								
			60	_								
			180	2.5.07	5,007							
			360	2.5627								
			1,440	1.1731	1 1.7791 2.4435 10.9325 8.0347 4.2311							

(continued)

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	osare to Racemic		VX- (ng/	·G				
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
			0	BDL	•	1	1	-			
			1	3.5232							
			10	3.1968	No samples received						
			20	3.5219							
			30	3.4886							
186			40	3.5935		110	samples receiv	veu			
			50	3.2308							
			60	3.0851							
			180	2.9320							
			360	2.4387							
		7.0	1,440	1.3646	1.6335 3.4071 26.7402 61.3619 5.35						
		7.0	0	BDL							
	30 Mar 2016		1	3.3067							
	30 Wai 2010		10	2.6189							
			20	2.7354							
			30	2.3520		N	o samples receiv	wed			
187			40	2.4331		110	samples recei	veu			
			50	2.3667							
			60	2.8746							
			180	2.0920	426						
			360	1.6426							
			1,440	0.9480	1.1248	1.7438	13.4322	45.1526	3.9924		
			0	118%*							
190		6.0	1	3.5114		No samples received					
170		0.0	10	3.2312					1		
21 3.1273 2.8458 9.9184 30.8							30.8527	54.4607	8.8111		

^{*}This sample was spiked with VX-G to determine percentage recovered from assay. BDL, below detection limit (<0.05 ng/g).

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time				K-G g/g)					
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain			
			0	BDL								
			1	2.7796								
			10	2.4796								
			20	2.4575								
			30	2.8872		N	o complec receiv	vad				
191			40	2.2457	No samples received							
			50	2.2163								
			60	2.2071								
		_	180	2.0524								
			360	1.7368								
	30 Mar 2016	6.0	1,440	0.9325	0.4247	0.8924	10.5911	30.5355	2.6535			
	30 Wiai 2010	0.0	0	BDL								
			1	3.0049								
			10	2.7175								
			20	3.0337								
			30	2.8918		N	o samples receiv	vad				
193			40	2.8594		IN	o samples fecely	eu				
			50	2.6646								
			60	2.6270								
			180	2.4253								
			360	1.9892								
			1,440	1.1554								

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	sure to Racellic		VX	K-G g/g)	,			
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
			0	BDL							
195			1	3.4135		N	o samples receiv	red			
			10	3.0525							
			20	2.7628							
			30	2.3390							
		6.0	40	2.4343	No samples received						
195			50	2.3165							
193			60	2.6175							
		-	180	2.1845							
			360	1.8254							
	30 Mar 2016		1,440	0.9291	1.4121	3.3995	23.8837	69.3294	4.2496		
	30 Wai 2010		0	BDL							
			1	4.5289							
			10	3.0838							
			20	3.3386							
			30	3.2996		N	o samples receiv	red			
196		5.0	40	2.8444		11	o samples receiv	cu			
			50	3.0163							
			60	2.9157							
			180	2.3040							
			360	2.4027		T	1				
I			1,440	1.0515	0.7072	1.3290	9.0184	14.0947	4.0463		

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	date to Ruceime		VX	X-G g/g)	,												
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain											
			0	BDL																
			1	2.9508																
			10	2.7473			No samples received													
			20	2.5805																
			30	2.5545		N														
197			40	2.3774		No samples received														
			50	2.5119																
			60	2.4025																
			180	2.0195																
			360	ISx																
	30 Mar 2016	5.0	1,440	1.0087	0.4971	1.0177	7.6585	20.6375	2.3159											
198	30 Wai 2010	3.0	0	BDL		N	o samples receiv	red												
			1	3.3165																
			10	3.0786																
				-	 	_	-	-							20	2.8316				
															- -				-	-
198	40 2.6164 No samples received																			
170			50	2.9105																
			60	2.8373																
			180	1.9559																
			360	1.6902			1	,												
			1,440	0.7815	815 0.5211 0.9505 8.2058 26.2812 3.5328															

(continued)

BDL, below detection limit (<0.05 ng/g). ISx, an insufficient amount of sample was received for assay.

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	VX-G (ng/g)							
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
		5.0	0	BDL							
			1	2.2156							
			10	2.1580							
			20	2.2090							
			30	2.3118	No samples received						
199			40	2.5210							
			50	2.0540							
			60	2.2927							
			180	2.0414							
			360	1.4989							
	30 Mar 2016		1,440	0.8579	0.4703	1.0219	.0219 7.5343 17.3610	17.3610	2.6870		
	30 Wiai 2010	4.1	0	BDL							
			1	2.4862							
			10	2.3898							
			20	2.3243							
			30	2.5103		No complex received					
200			40	2.3902	No samples received						
			50	2.3546							
			60	2.2639							
			180	2.0389							
			360	1.6037							
			1,440	1.0195	0.4193	4.6074	0.6096	8.3628	2.6363		

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	VX-G (ng/g)						
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain	
			0	BDL	-				•	
			1	2.5105						
			10	2.3213						
			20	2.3041						
			30	2.3237						
201			40	2.3818	No samples received					
			50	2.2290						
			60	2.2223						
			180	2.0126						
			360	1.5143						
	30 Mar 2016	4.1	1,440	0.9546	0.6681	0.8079	0.8079 2.8509 2.7339	2.7339	1.9049	
	30 Wai 2010	4.1	0	BDL						
			1	2.3394						
			10	2.3736						
			20	2.5656						
			30	2.9234		No complex received				
202			40	2.3611	No samples received					
			50	2.4053						
			60	2.1873						
			180	1.8610						
			360	1.7075						
			1,440	0.9754	0.5325	0.8058	8.2207	16.8110	3.3287	

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Lime	VX-G (ng/g)						
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver samples receive 6.3322	Kidney	Brain	
			0	BDL						
			1	2.2065						
			10	1.8914						
			20	2.1044						
203		30	1.8926	No someles resided						
203		4.1	40	1.9985	No samples received					
			50	1.8519						
			60	1.8539						
			180	1.6124						
			360	1.4551						
	30 Mar 2016		1,440	0.8176	0.4076	76 0.6700 6.3322 18.2030	18.2030	2.7183		
		3.7	0	BDL						
			10	2.2007						
			20	2.0657						
			30	1.9999						
204			40	1.8301	No samples received					
204			50	1.8139						
			60	1.4865						
			180	1.6931						
			360	1.5708						
			1,440	0.7391	0.2631	0.4417	2.7774	1.1685	1.5870	

Table 8. Results from VX-G Assays of Blood, Tissues, and Organs Following Intravenous Exposure to Racemic VX for Pharmacokinetic Studies (Continued)

Guinea Pig	Exposure	Dose	Sample Time	VX-G (ng/g)							
No.	Date	(µg/kg)	(min)	Whole Blood	Heart	Lung	Liver	Kidney	Brain		
				0	BDL			•			
			1	2.1550							
	10 2.4244 20 2.4614 30 2.3922 40 2.3494		10	2.4244							
			20	2.4614							
		No complex received									
206			40	2.3494 2.4830	eu						
			50	2.4830							
			60	2.1063							
			180	2.0363							
			360	1.6765							
	30 Mar 2016	3.7	1,440	1.0087	0.6088	1.3821	7.7758	18.9076	2.8424		
	30 Wiai 2010	3.7	0	BDL							
			1	3.2942							
			10	2.7747							
208			20	2.7421		N	No samples received				
			30	2.7562							
			40	2.4678							
			50	2.6672							
			60	2.4518							
208			180	2.3491		N	o samples receiv	ed			
200			360	2.3268							
			1,440	1.1310	0.4591	0.8752	6.9101	16.3958	3.2021		

BDL, below detection limit (<0.05 ng/g).




Figure 3. VX-G concentration in blood as a function of time for three different doses of P(+)-VX.

Table 9. Pharmacokinetic Parameters from Serial Samplings Following Intravenous Exposure to Various Doses of P(+)-VX, P(-)-VX, and a Racemic Mixture of VX

Stereoisomer	Dose (µg/kg)	k (min ⁻¹)	<i>t</i> _½ (min)
	245.0	0.0166	42
P(+)-VX	260.0	0.0099	70
	280.0	0.0354	20
	3.0	0.0011	635
D() WV	4.0	0.0010	694
P(-)-VX	4.2	0.0009	802
	5.0	0.0012	578
	3.7	0.0007	1,001
	4.1	0.0007	939
Racemic VX	4.8	0.0009	729
mixture	5.0	0.0010	720
	6.0	0.0009	803
	7.0	0.0008	825

Blank

LITERATURE CITED

- Bae, S.Y.; Winemiller, M.D. Separation of VX, RVX, and GB Enantiomers Using Liquid Chromatography—Time-of-Flight Mass Spectrometry; ECBC-TR-1341; U.S. Army Edgewood Chemical Biological Center: Aberdeen Proving Ground, MD, 2016; UNCLASSIFIED Report (AD1002738).
- Benschop, H.P.; de Jong, L.P. Nerve Agent Stereoisomers: Analysis, Isolation, and Toxicology. *Acc. Chem. Res.* **1988**, *21*, 368–374.
- Benschop, H.P.; de Jong, L.P. Toxicokinetics of Nerve Agents. In *Chemical Warfare Agents: Toxicity at Low Levels*; Somani, S.M.; Romano Jr., J.A., Eds.; CRC Press: Boca Raton, FL, 2001, pp 25–81.
- Maxwell, D.M.; Brecht, K.M.; O'Neill, B.L. The Effect of Carboxylesterase Inhibition on Interspecies Differences in Soman Toxicity. *Toxicol. Lett.* **1987**, *39* (1), 35–42.
- McGuire, J.M.; Demond, P.S.; Busch, M.W. *GC-MS/MS Analyses of Biological Samples in Support of Developmental Toxic Effects on Whole-Body Exposure of Rats to GB*; ECBC-TR-1286; U.S. Army Edgewood Chemical Biological Center: Aberdeen Proving Ground, MD, 2015; UNCLASSIFIED Report (ADA616556).
- McGuire, J.M.; Lumley, L.A.; Demond, P.S.; Busch, M.W.; Wright, L.K. *GC-MS/MS Analyses of Biological Samples in Support of Developmental Toxic Effects on Subcutaneous Exposure of Rats to GB*; ECBC-TR-1322; U.S. Army Edgewood Chemical Biological Center: Aberdeen Proving Ground, MD, 2015; UNCLASSIFIED Report (ADA624375).
- Nordgren, I.; Lundgren, G.; Puu, G.; Holmstedt, B. Stereoselectivity of Enzymes Involved in Toxicity and Detoxification of Soman. *Arch. Toxicol.* **1984,** *55* (1), 70–75.
- Ordentlich, A.; Barak, D.; Sod-Moriah, G.; Kaplan, D.; Mizrahi, D.; Segall, Y.; Kronman, C.; Karton, Y.; Lazar, A.; Marcus, D.; Velan, B.; Shafferman, A. Stereoselectivity toward VX Is Determined by Interactions with Residues of the Acyl Pocket as well as of the Peripheral Anionic Site of AChE. *Biochemistry* **2004**, *43* (35), 11255–11265.
- Shih, T.M.; McDonough, J.H. Efficacy of Biperiden and Atropine as Anticonvulsant Treatment for Organophosphorus Nerve Agent Intoxication. *Arch. Toxicol.* **2000**, *74*, 165–172.
- Wright L.K.; Forster J.S.; Moretz R.W.; Gaviola B.I.; Renner J.A.; Kristovich R.L. *Median Lethal Doses Associated with Intravenous Exposure to the Optically Pure Enantiomers of VX in Guinea Pigs*; ECBC-TR-1437; U.S. Army Edgewood Chemical Biological Center: Aberdeen Proving Ground, MD, 2017; UNCLASSIFIED Report.

Blank

ACRONYMS AND ABBREVIATIONS

BDL below detection limit

CAS Chemical Abstracts Service

CE collision energy CI chemical ionization

ECBC U.S. Army Edgewood Chemical Biological Center

GC gas chromatography

IACUC Institutional Animal Care and Use Committee

IPA 2-propanol

IS internal standard

ISxinsufficient amount of samplekelimination rate constantKFpotassium fluorideLD50median lethal dose

MRM multiple reaction monitoring

MS mass spectrometry

MS/MS tandem mass spectrometry m/z mass-to-charge ratio NB laboratory notebook

NMR nuclear magnetic resonance R^2 coefficient of determination SPE solid-phase extraction

t√₂ half-life

VX O-ethyl S-(2-diisopropylaminoethyl) methylphosphonothioate

VX-G O-ethyl methylphosphonofluoridate, EA 1207

DISTRIBUTION LIST

The following individuals and organizations were provided with one Adobe portable document format (pdf) electronic version of this report:

U.S. Army Edgewood Chemical Biological Center (ECBC) RDCB-DRT-T ATTN: McGuire, J.

Evans, R.

Defense Threat Reduction Agency

J9-CBS

ATTN: Graziano, A.

Office of the Chief Counsel

AMSRD-CC

ATTN: Upchurch, V.

G-3 History Office U.S. Army RDECOM

ATTN: Smart, J.

Department of Homeland Security

RDCB-PI-CSAC ATTN: Negron, A. DHS-S&T-RDP-CSAC

ATTN: Strang, P.

ECBC Rock Island

RDCB-DES ATTN: Lee, K. RDCB-DEM

ATTN: Grodecki, J.

ECBC Technical Library

RDCB-DRB-BL ATTN: Foppiano, S. Stein, J.

Defense Technical Information Center

ATTN: DTIC OA

